


# NSV Instruments

## Li Brand™ No-Scalpel Vasectomy Instruments

- Used and proven in every international published study on NSV
- Made by a German company to AMI's exacting standards
- Come with a No-Risk, limited 5-year warranty

Li Brand™ NSV instruments are the Original NSV instruments with the superior cantilevered design proven in literally millions of cases. Independent quality evaluation by Indeps, Inc. has shown that Li Brand™ instruments outperform and last longer than others. It's no wonder that so many top urologists and major organizations choose Li Brand™ and Advanced Meditech as their supplier. Get the very best. Get Li Brand™ instruments.



### Did You Know?

The cantilever on Li Brand™ ringed forceps is the critical factor in NSV's documented 10-fold reduction in complications, especially hematoma. As the Li Brand™ ringed forceps is locked around the vas extracutaneously, its tips spring slightly apart (Figure 1). As a result, the forceps holds firmly without pinching or crushing the delicate tissues around the vas. (Figure 2)

Figure 1

Figure 2

# NSV Education

## QuickStart™ NSV Training Exclusives

Advanced Meditech is the exclusive distributor of *No-Scalpel Vasectomy: The Goldstein Video* (VED-2) with Marc Goldstein, M.D. of The New York Presbyterian Hospital – Weill Medical College of Cornell University. With superb surgical footage, Dr. Goldstein provides step-by-step instruction to help you master NSV.

With its exclusive QuickStart™ program, AMI also loans no-scalpel vasectomy training instruments to residency training programs at little or no cost. For more information, send an inquiry to [info@ameditech.com](mailto:info@ameditech.com).

Advanced Meditech also provides many other training and educational products for physicians and patients, including the *NSV Manual for Surgeons* (VED-1), patient education video (VED-3), scrotal training model (VED-9, shown at right), and our best-selling Personal-Physician™ NSV Patient Education Diagrams and Brochures VED-6 and VED-7).



### Did You Know?

Hands-on training is the fastest, most efficient way to achieve mastery of NSV. With hands-on training, practitioners report that they can reduce procedural time by up to 50%, while reducing or eliminating most operative pain and bleeding.

