


NSV Counseling

Well-Presented Information Means Higher Patient Satisfaction and Increased Referrals

- VED-3 No-Scalpel Vasectomy Patient Education Video 10 minutes, English, VHS. Discusses making the decision to have a vasectomy. Explains vasectomy and NSV. Great for pre-consults.
- VED-5 No-Scalpel Vasectomy: Good News for Men 8pp., Brochure answering commonly asked questions about no-scalpel vasectomy.
- VED-6 PersonalPhysician™ NSV Counseling Brochures 6pp., 4/c tri-fold, CUSTOMIZED with individual physician and practice information, even logos. Excellent referral generator. Also in Spanish VED-6S.
- VED-7 Full Color NSV Male Reproductive Anatomy Diagram 2-sided 8.5"x11" 4/c laminated illustration, male anatomy on one side; anatomy and NSV on the other. Designed for use with VED-6 brochures or as a stand alone counseling tool.



Clockwise from top: VED-6 PersonalPhysician Custom Brochures, VED-7 Diagram, VED-5 Brochures, and VED-3 Patient Education Video.

Did You Know?

No-Scalpel Vasectomy (NSV) is rapidly displacing conventional incisional vasectomy in the United States. Despite being introduced to the U.S. less than 15 years ago, approximately 170,000 procedures, about 1/3 of all vasectomies in the U.S., are now performed annually using the no-scalpel technique.

NSV Accessories

Give You Easier Procedures and Longer Instrument Life

- AC-7540 Sterilizable Instrument Case Autoclavable, 7.5"x4"x0.75", with silicone insert
- VE-4 Iris Scissors German-made, 4.5" long, straight.
- VE-5 Surgical Instrument Lubricant Anticorrosive, steam compatible. Helps keep your Li Brand™ instruments looking like new
- VE-11 Handzfree™ Anesthetic Bottle Holder Holds 4 x 50ml bottles on wall or in cabinet; allows drawing of anesthetic without an assistant



Clockwise from top: VE-11 Handzfree™ Anesthetic Bottle Holder, VE-4 Iris Scissors, VE-5 Instrument Lubricant, AC-7540 Case with NSV Instruments

Did You Know?

Proper instrument handling and care can extend the life of your NSV instruments by up to 500%. NSV dissecting forceps have tips as delicate as those of microsurgical instruments. Unlike other sharp instruments, NSV dissecting forceps *should never be placed tip down* on metal trays or other hard surfaces.